



ULTRADENT

Thermo Clone VPS Impression Material

Ultradent's Thermo Clone VPS impression material represents an exciting advancement in prosthodontics and restorative dentistry. To the dental professional taking an impression while trying to minimize patient discomfort, bubble gum-scented Thermo Clone is the brand of vinyl polysiloxane impression material that makes impression-taking easier because of its consistent setting times, thermal set-curing that activates when placed in the patient's mouth and the most effective techniques to getting accurate impressions.

Thermo Clone VPS impression material is specially engineered to be reliable and consistent. Its extremely temperature-sensitive nature allows for long working time, a short setting time and minimal time in the critical zone where distortions occur. The product allows for up to a minute and a half of working time (compared to the leading competitor's 35 seconds of working time) to reduce clinician error, and it sets more than 2 minutes faster in the patient's mouth (at 2.25 minutes, compared to an average of 4 minutes with the leading competitor), leading to fewer distortions because many patients inadvertently move and grow uncomfortable if the impression material is in their mouth for too long.

Thermo Clone VPS impression material is also hydrophilic, ensuring precision in every impression by displacing moisture so

the material can reach every part of every tooth. The material's high tear strength allows clinicians to remove it confidently from a patient's mouth without being worried about breaks or tears, while also demonstrating excellent dimensional stability, maintaining its shape during transfer to the dental lab.

Thermo Clone VPS impression material works with all automatic mixing machines and comes with Ultradent's unique dynamic mixing tip.

Time—literally every single minute—is incredibly valuable to clinicians. When dentists use an impression material that features consistent setting times and a thermal-accelerated set, such as Thermo Clone material, they put valuable, hard-earned money back in their pockets while getting the right impression the very first time. ■

For more information, call 800-552-5512
or visit ultradent.com.

ULTRADENT
PRODUCTS, INC.