

GARRISON DENTAL

COMPOSI-TIGHT 3D FUSION SECTIONAL MATRIX

Garrison's Composi-Tight 3D Fusion Sectional Matrix is engineered to successfully complete a wider variety of Class II composite restorations. This comprehensive kit includes three varieties of separator rings, five sizes of matrix bands, four sizes of interproximal wedges and a ring placement forceps. Composi-Tight 3D Fusion is indicated for use in posterior Class II MO/DO, MOD and wide preparation restorations.

Garrison
Dental Solutions

For more information, visit garrisondental.com.



HIGHLIGHTS

- The three-ring system consists of **short, tall and wide preparation separator rings**. All three incorporate Garrison's Soft-Face silicone tips for flash reduction and Ultra-Grip retention extensions that eliminate ring "spring-off."
- Five sizes of ultrathin, dead-soft nonstick sectional matrix bands have a Grab-Tab on the occlusal edge to **simplify band placement and removal**. The bands are coupled with four sizes of Composi-Tight 3D Fusion Ultra-Adaptive wedges.