



## **CAO GROUP**

## MONET LASER CURING LIGHT

aser curing has some major advantages over LED curing. The Monet Laser Curing Light produces a deeper cure than LED or any other light source during the same exposure time. The collimated beam and consistent power of the Monet make for superior bond strength and faster, deeper and more reliable curing results. With the proper application of the Monet laser, clinicians don't have to worry about weak bonds or soft cures. Each click of the Monet is a one-second 2.5 mm cure; bulk cures are possible in just three seconds!

The Monet has been tested by multiple independent sources with impressive results. Gordon J. Christensen's *Clinicians Report* testing revealed "high irradiance output produced faster and deeper polymerization than other light sources." *Clinicians Report* also tested the laser curing for heat and found zero heat issues when used properly. The Monet also comes with unique apertures to reduce the size or even reduce the power by 50% for even more control when precision is needed.

Sixty to 80 percent of dental procedures involve composite curing. By eliminating soft cures and cutting chair time, the Monet laser curing light improves the patient experience.

## **HIGHLIGHTS**

- A 13.2-mm-wide lens, a rotating 90-degree head angled for access to all tooth surfaces, and a charging station with an extra rechargeable battery.
- Comes with two safety goggles, glasses insert, a protective shield and paddle, four apertures and barrier sleeves.
- Consistent power to produce a deep and lasting cure in just one second with minimal shrinkage and no postoperative sensitivity.



For more information, visit

caogroup.com/products/monet-curing-laser.