

RESTORATIVE DENTISTRY

GINGIVAL RETRACTION MATERIALS

Ultradent Products: Ultrapak Retraction Cord

Ultrapak cord is made of 100% cotton, knitted into thousands of tiny loops to form long, interlocking chains. This unique knitted design exerts a gentle, continuous outward force following placement as the knitted loops seek to open. Optimal tissue displacement occurs in three to eight minutes.

**MULTIYEAR
WINNER**
since 2003



TG'S TAKE: *"Ultradent has a thick catalog filled with materials and equipment, and some of its oldest members of that catalog are still winning awards. The method for making Ultrapak retraction cord is so special that if you take a tour of the Ultradent facility, this room is off-limits! Even in the era of digital dentistry, tissue must be managed for a crisp, clear view of the finish line. Retraction cord remains as one of the best options for achieving this goal."*

130 OF 422 VOTES

Information: 888-230-1420 or ultradent.com

FREE FACTS: circle 91 on card.

NEW PRODUCTS Introduced between May 1, 2022 - April 30, 2023

NEW PRODUCTS: EQUIPMENT

Ultradent Products: Valo X Curing Light

46 OF 465 VOTES

Information: 888-230-1420
or ultradent.com

FREE FACTS: circle 92 on card.

NEW PRODUCTS: MATERIALS

Tokuyama Dental America: Omnichroma Flow Bulk

70 OF 471 VOTES

Information: 877-378-3548
or omnichromaflowbulk.com

FREE FACTS: circle 93 on card.